DESCRIZIONE:

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CE 0476

DESCRIPTION:

Drill bits are instruments that have to be used by surgeons and medical staff only. Drill bits are in stainless steel and are available in various length and diameter. Different type of Titanium (AISI 316 L, AISI 420 etc.) can be used for drill bits production.

THE MENTIONED MATERIALS IS PROVIDED DECONTAMINATED, UNSTERILIZED;

ACCESSORIES:

Drill bits have to be normally used in conjunction with electric equipments, not produced by Cizeta Surgical, and that are already in use in operating rooms. Using the conjunction of electric equipments with Cizeta Surgical drill bits guarantees in any case to comply with 93/42/CCE Directive.

INDICATIONS:

Drill bits are intended to be used for performing bone osteotomy or bone re-modelling; are also used to make holes for screws application.

CONTRAINDICATIONS:

1. Ascertained allergies to stainless steel.
2. Patient's special conditions: senility, alcoholism and infections

PRECAUTIONS:

Surgeon must not start clinical use of drill bits without having previously read the instructions manual; in addition, these devices are to be utilized exclusively in specialized environment (ambulatory or operating room). Medical staff must check the package integrity, before using the device and at the end of the sterilization process too. Cizeta Surgical it is not responsible for improper use of these devices.

Do not use if packaging is broken or damaged.

MAINTENANCE:

Devices have to be properly disinfected, cleaned and sterilized prior to be used on the patient.

Place instruments in a disinfectant solution containing special additives to prevent corrosion. Do not use harsh chemicals (eg. Hydrochloric acid, hydrogen peroxide), they may corrode instruments.

Clean the instruments in an ultrasonic cleaner or distilled water.

Before sterilization, instruments have to be stored in an adequate container.

After sterilization, instruments have to be kept in close recipients.

WARNINGS:

Ensuring a supply of water cooling with saline irrigation to prevent damage to the bone. The cooling liquid has to be properly distributed to the entire working surface of drill bit.

Devices with damaged cutting efficiency, or bent, or with an eccentric rotation should be immediately discarded and not used again as they may cause severe injury to both the user and patient.

Working pressure should be between 0,3 and 2,0 N/mm². Avoid excessive pressure.

Do not wrap or leverage on the device during processing.

Respect the maximum speed of rotation of 500 – 800 rev/min.

In case of using the optional extension, ensure drill bit by tightening the locking screw.

STERILIZATION:

Great attention is paid to the production of stainless steel drill bits in order to guaranty reliability and quality of operating outcome to surgeons. Medical staff too must contribute to achieve the expected outcome, paying due attention in handling and utilizing the devices. This set of instruments do not require special indications in sterilization method; Cizeta Surgical recommends autoclave sterilization, in a vacuum, at 134°C (273° F) temperature for a period of 10 minutes.

The chosen sterilization process, in any case, must comply with the EN556 norm which states that the theoretical probability of vital micro-organism presence on the device must be less or equal to 1 over 10⁶ (S.A.L. = 10⁻⁶, where S.A.L. stands for "Sterility Assurance Level"). It is the user's responsibility to ensure the proper sterilization process and to verify the barrenness of all devices and during the entire phases.

These device are assembled in clean room and supplied, **unsterilized**, in medical packs, one side made in medical grade paper and the other side in polypropylene/polyester transparent plastic. The packing bears two chemical indicators:

- for steam ⇒ dye toning to brown color change
- for ethylene oxide ⇒ dye toning to yellow color change

USAGE INSTRUCTION:

Sterilization is to be done following above instructions

Drill bits are to be handled exclusively in suitable environment (specialized ambulatory or operating room) and with due care (not to be used without sterilized gloves).

Program the speed of rotation more suitable to the type of surgical procedure to be performed.

Check the perfect correspondence between the drill bit with the corresponding electric equipment.

In case of re-use accurate cleaning, before sterilizing, is requested.

CIZETA SURGICAL S.r.l.
Via Caselle 76
40068 San Lazzaro di Savena (BO)
Tel: +39 051.462000

WARNING
REFER TO
INSTRUCTIONS
TO USE

AVOID
MOISTURE

CE 0476