DESCRIZIONE:

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DESCRIPTION:

Screwdrivers have a Teflon handle and a stainless steel blade and, as an optional, a centring device, which differs in relation to screw head type. Screwdrivers blades are carefully checked in order to ensure adequate screw retaining.

THE MENTIONED MATERIALS IS PROVIDED DECONTAMINATED, UNSTERILIZED.
THE MENTIONED MATERIAL CAN BE RE-USED AFTER PROPER MAINTENANCE (see specific chapter).

INDICATION:

Screwdrivers are constituted by a stainless steel rod for different series of screws and pins, and by a Teflon handle. This set of instruments is compact and easy to use; it can be easily stored, at hospital personnel or surgery staff care; this set of instruments can be reused after an adequate washing and relative sterilization.

CONTRAINDICATIONS:

1. Ascertained allergies to stainless steel
2. Patient's special conditions: senility, alcoholism and infections

PRECAUTIONS:

Surgeon must not start clinical use of screwdrivers without having previously read the instructions manual.

MAINTENANCE:

Devices have to be properly disinfected, cleaned and sterilized prior to be used on the patient. Place instruments in a disinfectant solution containing special additives to prevent corrosion. Do not use harsh chemicals (eg. Hydrochloric acid, hydrogen peroxide), they may corrode instruments. Clean the instruments in an ultrasonic cleaner or distilled water. Before sterilization, instruments have to be stored in an adequate container. After sterilization, instruments have to be kept in close recipients.

STERILIZATION:

Cizeta Surgical recommends autoclave sterilization, in a vacuum, at 134°C (273° F) temperature for a period of 10 minutes.

The chosen sterilization process, in any case, must comply with the EN556 norm which states that the theoretical probability of vital micro-organism presence on the device must be less or equal to 1 over 10^6 (S.A.L. = 10^{-6} , where S.A.L. stands for "Sterility Assurance Level"). It is the user's responsibility to ensure the proper sterilization process and to verify the barrenness of all devices and during the entire phases.

These device are assembled in clean room and supplied, **unsterilized**, in medical packs, one side made in medical grade paper and the other side in polypropylene/polyester transparent plastic. The packing bears two chemical indicators:

- for steam \Rightarrow dye toning to brown color change
- for ethylene oxide \Rightarrow dye toning to yellow color change

USAGE INSTRUCTION:

Sterilization is to be done following above instructions. These instruments are to be utilized exclusively in specialized environment (ambulatory or operating room).

Insert the screwdriver blade into the handle, by verifying that blade is very well secured.

Use the centring device (optional) to centre the head screw: place the centring device on screw head while is in the Teflon container; push the screwdriver blade up to hear the retaining click. Start then to use the device; at the end, the centring device will automatically disengage.