DESCRIZIONE:

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CE 0476

DESCRIPTION:

Titanium plates are implantable surgical devices, utilized in cranio-maxillo-facial surgery. Their various different shapes and extreme adaptability allow to comply with the anatomic shapes of cranial-facial bones as much as possible. Such devices are disposable and not reusable. The hereby reported instructions are valid for all plates, even in different forms and thickness, which are included in the following modular systems:

MODUL SYSTEM 1.0
MODUL SYSTEM 1.2
MODUL SYSTEM 1.6
MODUL SYSTEM 2.0
MODUL SYSTEM 2.3
MODUL SYSTEM 2.7

ABOVE MENTIONED MATERIALS ARE PROVIDED DECONTAMINATED, UNSTERILIZED, DISPOSABLE; NOT REUSABLE

INSTRUCTIONS:

- Synthesis of Cranial mandibular fractures
- Correction of bone malformation
- Mandibular reconstruction

WARNINGS:

- Ascertained allergies to titanium
- Patient's special conditions: senility, alcoholism and infections
- Device reuse

The devices are disposable and cannot be re-sterilized.
Potential risk in case of re-use: INFECTION, DEFORMATION.

PRECAUTIONS:

Surgeon must not start clinical use of plates without having previously read the instruction manual; in addition, plates are to be used in specialized environment (ambulatory or operating room). Medical staff must check the package integrity, before using the device and at the end of the sterilization process

ATTENTION:

Titanium plates are furnished decontaminated but not sterilized. Therefore, before utilization, plates must be sterilized and handled with due precaution to avoid infections.

STERILIZATION:

Great attention is paid to the production of plates in order to guarantee reliability and quality of operating outcome to surgeons. Medical staff too must contribute to achieve the expected outcome, paying due attention in handling and utilizing the devices. Titanium plates do not require special indications in sterilization method; Cizeta Surgical recommends autoclave sterilization, in a vacuum, at 134°C (273° F.) for a period of 10 minutes.

The chosen sterilization process, in any case, must comply with the EN556 norm which states that the theoretical probability of vital micro-organism presence on the device must be less or equal to 1 over 10⁶ (S.A.L. = 10⁻⁶, where S.A.L. stands for "Sterility Assurance Level"). It is the user's responsibility to ensure the proper sterilization process and to verify the barrenness of all devices and during the entire phases.

These device are assembled in clean room and supplied, **unsterilized**, in medical packs, one side made in medical grade paper and the other side in polypropylene/polyester transparent plastic. The packing bears two chemical indicators:

- for steam ⇒ dye toning to brown color
- for ethylene oxide ⇒ dye toning to yellow color

USAGE INSTRUCTION:

1. The plate's sterilization is to be done following above instructions.
2. Plates are to be handled exclusively in suitable environment (specialized ambulatory or operating room) and with due care (not to be used without sterilized gloves).
3. Plates are to be applied and adapted as the situation may require
4. Device is to be utilized with the screws outfit

For device application, Cizeta Surgical "Modul System" proper screws are required.

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CE 0476

DESCRIPTION:

Titanium screws and pins are implantable surgical devices, utilized in cranio-maxillo-facial surgery to fix the relevant titanium plates. Titanium is the ideal implantable material, well known for its excellent biocompatibility, resistance to corrosion and no-toxicity in biological environment.

Screws are made in various shapes and lengths: the screws common diameter is: 1.0, 1.2, 1.6, 1.8, 2.0, 2.3, 2.7, 3.0 ; length varies from 2 mm. to 30 mm.

Pins are available in 1,0mm diameter and lengths of 3,0 or 5,0 mm.

THE MENTIONED MATERIALS IS PROVIDED DECONTAMINATED, UNSTERILIZED, DISPOSABLE; NOT REUSABLE.

INDICATION:

Exclusively for the synthesis of relevant plates, or bone and membrane fragments and bone distractors.

CONTRAINDICATIONS:

1. Ascertained allergies to titanium
2. Patient's special conditions: senility, alcoholism and infections
3. Device reuse

The device is disposable and its reuse involves the following major issues:

- INFECTION: even subsequent manipulation does not guarantee a microbial conform to the type of sterilization. Possible pyrogenicity.
- STRAIN: the device could suffer in the subsequent manipulation of small deformations that would make it unsuitable for use.
- Re-sterilization could lead to serious infections.

PRECAUTIONS:

Surgeon must not start clinical use of screws and or pins without having previously read the instructions manual; in addition, screws, pins and relevant instruments are to be utilized exclusively in specialized environment (ambulatory or operating room). Medical staff must check the package integrity, before using the device and at the end of the sterilization process too.

ATTENTION:

SCREWS AND PINS ARE PROVIDED DECONTAMINATED UNSTERILIZED.

Screws and pins must be sterilized and handled with due precautions in order to avoid infections.

STERILIZATION:

Great attention is paid to the production of titanium screws and pins in order to guaranty reliability and quality of operating outcome to surgeons. Medical staff too must contribute to achieve the expected outcome, paying due attention in handling and utilizing the devices. These devices do not require special indications in sterilization method; Cizeta Surgical recommends autoclave sterilization, in a vacuum, at 134°C (273° F) temperature for a period of 10 minutes.

The device can be sterilized with ethylene oxide also. The sterilization process must be validated in accordance with the UNI EN ISO 11135-1:2008. Warning: the residue of ethylene oxide, if not intended, can cause serious infections

The chosen sterilization process, in any case, must comply with the EN556 norm which states that the theoretical probability of vital micro-organism presence on the device must be less or equal to 1 over 10^6 (S.A.L. = 10^{-6} , where S.A.L. stands for "Sterility Assurance Level"). It is the user's responsibility to ensure the proper sterilization process and to verify the barrenness of all devices and during the entire phases.

These devices are assembled in clean room and supplied, **unsterilized**, in medical packs, one side made in medical grade paper and the other side in polypropylene/polyester transparent plastic. The packing bears two chemical indicators:

- for steam \Rightarrow dye toning to brown color change
- for ethylene oxide \Rightarrow dye toning to yellow color change

USAGE INSTRUCTIONS:

SCREWS AND PINS: Sterilization is to be done following above instructions.

1. Screws and pins are to be handled exclusively in suitable environment (specialized ambulatory or operating room) and with due care (not to be used without sterilized gloves).
2. Screws are to be applied after having previously drilled the bone with the relevant drill bits accordingly to the screw length; be sure also to apply less screwing torsional force than traditional single thread screws; afterwards, the proper plates can be applied.

CIZETA SURGICAL S.r.l.

Via Caselle 76

40068 San Lazzaro di Savena (BO)

Tel: +39 051.462000

SINGLE USE

WARNING
REFER TO
INSTRUCTIONS
TO USE

AVOID
MOISTURE

CE 0476