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TEME L'UMIDITA'

MEDICO

INSTRUCTIONS FOR USE TITANIUM SCREWS AND NAILS

Rev.05 del 08/23 Pag. 1/3

DESCRIZIONE:

Titanium screws and pinls are implantable surgical devices used in cranio- maxillofacial area surgery and to attach corresponding titanium plates or bone fragments to each other or prostheses in general. Titanium is the ideal implantable material, known for its excellent biocompatibility, corrosion resistance, and nontoxicity in a biological environment. Material: Screws and nails are made of Titanium Grade 5 Ti6Al4V - ASTM F136 and ISO 5832-3:

CHEMICAL COMPOSITION OF THE MATERIAL:

Element	Maximum composition
	<u>limits</u>
	Mass fraction %
Aluminum	5.5 to 6.75
<u>Vanadium</u>	3.5 to 4.5
Iron	Maximum 0.3
<u>Oxygen</u>	Maximum 0.2
<u>Carbon</u>	Maximum 0.08
Nitrogen	Maximum 0.05
Hydrogen	Maximum 0.015
<u>Titanium</u>	Balance

There are screws with different head shapes (square, cross, star or torx engagement); the types used are with diameters of 1.0, 1.2, 1.4, 1.6, 1.8, 2.0, 2.3, 2.7, 3.0: lengths range from 2 mm to 30 mm. Pins and screws nail 1.0 mm in diameter and 3.0 or 5.0 mm lengths.

THE ABOVE MATERIAL IS SUPPLIED CLEAN, NONSTERILE, DISPOSABLE; IT IS NOT REUSABLE.

INTENDED USE:

Exclusively for attachment of related plates, membranes or MESH, bone fragments to each other, distractors, wires for intermaxillary locking or orthodontic anchorage.

USE OF THE PRODUCT IS STRICTLY RESERVED FOR SURGEONS AND/OR BY SURGICAL STAFF IN MAXILLOFACIAL FRACTURE SYNTHESIS AND MANDIBULAR RECONSTRUCTIONS.

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- Established allergies to titanium.
- Special patient conditions such as senility, alcoholism, and infection.
- Reuse of the device.

The device is disposable, and its reuse leads to the following serious problems:

- INFECTION: Even subsequent handling does not guarantee a bacterial load in accordance with the type of sterilization. Possible pyrogenicity and failure to sterilize resulting in infection.
- DEFORMATION: The device may undergo in subsequent handling small deformations that would make it unfit for use.

INSTRUCTIONS FOR USE TITANIUM SCREWS AND NAILS

Rev.05 del 08/23 Pag. 2/3

PRECAUTIONS:

The surgeon should initiate clinical use of the screws and pins only after having first read the instructions for use; in addition, he or she should use these devices and related instrumentation in a specialized setting (outpatient clinic or operating room). Medical personnel should check the integrity of the packaging before using the product and upon completion of sterilization.

WARNING:

SCREWS AND PIN S ARE SUPPLIED CLEAN BUT NOT STERILIZED.

Screws and pinls, before their use, should be sterilized and handled with proper precautions to avoid contamination. **The screws are made using the double-principle fast-forward technique**. This feature requires half the torque force for screwing the screw onto the bone. You must report any serious incidents that have occurred in connection with the medical device we supply to the manufacturer and to the competent authority of the member state in which you are located.

STERILIZATION:

Titanium screws and pinls are manufactured with great care to ensure the surgeon's safety and quality of the surgical result. Medical personnel must also contribute to the achievement of the result by taking due care in handling and using the devices. There are no special restrictions in the sterilization method for these devices; Cizeta recommends autoclaving at a temperature of 134°C for a duration of 10 minutes.

Beware if sterilized with ethylene oxide: the residue, if not within legal limits, can be harmful. Remember that the sterilization process chosen must still be able to ensure compliance with EN556, which states that the theoretical probability of viable microorganisms being present on the device must be less than or equal to 1 in 10_6 (S.A.L. = 10_{-6} , where S.A.L. means Sterility Assurance Level). It is the responsibility of the user to validate all sterilization steps and verify sterility for all products and in all steps used. It should be noted that once the product has been sterilized it cannot be re-sterilized. The present devices, supplied **non-sterile**, are packaged in bags having one side made of medical grade paper and one side made of transparent polypropylene/polyester plastic. There are two chemical indicators on the package including one for one for vapor \Rightarrow brown-colored viraggio and one for ETO yellow-colored viraggio.

INSTRUCTIONS FOR USE:

SCREWS AND PINS: Carry out sterilization of screws and nails according to the instructions above.

Handle screws and nails only in a suitable environment (specialized clinic or operating room) with proper precautions (do not touch them without sterile gloves). Proceed with the application of the screws after drilling the entry hole in the bone with the relevant drill, taking care to choose the appropriate length of the drill in relation to the length of the chosen screw; apply less twisting force in the assembly of the screw clamping on plate than in traditional single-principle screws. Then use the corresponding plates if the need arises.

CONSERVATION

Store in a way that prevents damage to the package. Moisture tolerant. Store at room temperature in a dry place. Not reusable. The expiration date shown refers to the product in unopened package. Do not use the device if the packaging is damaged.

DISPOSAL

The devices are used in a hospital setting therefore they must be disposed of according to current waste disposal laws and enforced by the hospital facility. In particular, devices that are used and therefore may be polluted and biologically contaminated must be disposed of as special waste. Since they are intended to be used at hospital facilities that routinely implement such special collection, it is not necessary to provide such indications in labeling.

The safety and clinical performance summary is available from EUDAMED (https://ec.europa.eu/tools/eudamed) and until that platform is fully operational can be requested from CIZETA SURGICAL SRL by sending an email to mmarfori@cizetasurgical.it

INSTRUCTIONS FOR USE TITANIUM SCREWS AND NAILS

Rev.05 del 08/23 Pag. 3/3

CIZETA SURGICAL S.r.l.

Via Caselle 76 40068 San Lazzaro di Savena (BO)

MEDICAL DEVICE

READ THE INSTRUCTIO NS FOR USE