

**Distrattori Alveolari**  
**Distrattori Palatali**  
*Alveolar distractors*  
*Palatal distractors*

## Descrizione

Apparecchio per la distrazione verticale della sede alveolare:

Per ottenere la corretta verticalità di distrazione del tassello osseo occorre spessorare l'ancoraggio del corpo mobile con placche Cizeta Surgical della serie Modul System 1,6.

- Disponibile in tre versioni con ampiezza di distrazione fino a 10 mm, 15 mm e 20 mm.
- Poco invasivo per il ridotto spessore del dispositivo, 2 mm.
- Da impiantare con viti di Ø 1,6 Cizeta Surgical e relativo strumentario.

## Description

*Distraction device to uplift the alveolar site:*

*In order to obtain the proper distraction verticality of the bone segment, the anchorage of the mobile body must be thickened with Cizeta Surgical plates plates from the Modul System 1,6 set.*

- *Available in three patterns with distraction range up to 10 mm, 15 mm and 20 mm.*
- *Low profile due to device's 2 mm thickness.*
- *To be implanted with Ø 1,6 Cizeta Surgical screws and related instruments.*

