

Distrattori Alveolari
Distrattori Palatali
Alveolar distractors
Palatal distractors

Descrizione

Apparecchio per la distrazione del palato:

La vite di distrazione, con arresto a fine distrazione, non è scomponibile e, opportunamente ancorata ad un dente o a una vite, non può essere ingerita dal paziente. La placca di alloggiamento si fissa con un pin e con una vite dalla parte vestibolare.

- Disponibile in due possibilità di distrazione:
da 18 mm a 30 mm
da 28 mm a 45 mm
- Da impiantare con viti di Ø 2,0 Cizeta Surgical e relativo strumentario.

Description

Distraction device of the palate:

The distraction screw is made in such a way that at its maximum opening, it cannot be disassembled; appropriately anchored to a tooth or a screw, it cannot be ingested by the patient. The fixing plate is fastened to vestibular part by a pin and a screw.

- *Available in two options of distraction:
from 18 mm to 30 mm
from 28 mm to 45 mm*
- *To be implanted with Ø 2,0 Cizeta Surgical screws and related instruments.*

