



**Distrattori**  
**Cranio - maxillo - facciali**  
*Cranio - maxillofacial distractors*

## Descrizione

Apparecchio per la distrazione orbitale con asta di attivazione extra tessuti molli.

- Applicabile alla distanza richiesta dall'orbita del paziente;
- Agisce con una forchetta per attivare la distrazione senza la necessità di fissare la parte da distrarre con viti;

## Vantaggi

- A distrazione ultimata si rimuovono le viti di fissaggio sul cranio e la forchetta senza intervenire chirurgicamente sull'orbita.
- Si fissa con viti di Ø 2,0 Cizeta Surgical e relativo strumentario.

## Description

Device for orbital distraction with extra soft tissue activation bar.

- Can be positioned at the required distance from the orbit;
- It works with a fork to activate the distraction avoiding the need of screwing the distraction area;

## Benefits

- At the end of distraction period, both fixing cranial screws and orbital fork can be removed without surgery in orbital area.
- To be implanted with 2,0 Cizeta Surgical screws and related instruments.

