

Contenitori Sterilizzabili

Sterilizable containers

Descrizione

I contenitori appositamente realizzati per accomodare lo strumentario chirurgico Cizeta Surgical sono realizzati in Alluminio anodizzato e/o in acciaio inox con leve di chiusura in acciaio inossidabile che consentono una apertura aseptica del coperchio.

Vantaggi

- Optional, per una facile identificazione del container, sono disponibili su di un lato due porta targhette in cui è possibile inserire delle targhette in alluminio;
- La valvola infine è studiata appositamente affinché rimanga chiusa in condizioni di utilizzo normali, mentre si apre in condizioni di sovra o sotto pressione.

Description

The containers are designed for properly accommodate Cizeta Surgical instruments; they are made in Anodized Aluminum and/ or stainless steel with stainless steel closing levers that guarantee a contamination-free removal of the lid.

Benefits

- To ensure easy container identification, there is the possibility (optional) of inserting aluminum label into two label holders;
- Finally the upper valve is a device designed to remain closed in normal condition, while in the event of excess or low pressure, it opens.

05CO-SE Contenitore sterilizzabile 28 x 28 cm
05CO-SB Sterilizable containers 28 x 28 cm



05CO-SE2 Contenitore sterilizzabile 30 x 13
05CO-SB2 Sterilizable containers 30 x 13

